

***Harvesting Protein Lysates from 3-D Acinar Cultures  
(for quantification of non-phosphorylated, intracellular, stable proteins)***

- 1) Aspirate the media from the wells for harvesting.
- 2) Add high concentration trypsin (0.25%) for washing and remove immediately. For a single well of an 8-well chamber slide, use 250  $\mu$ l of trypsin for washing. For cultures grown in larger wells, such as 24-well or 12-well plates, use 500  $\mu$ l or 1 ml, respectively for washing.
- 3) Add half of the volume of high concentration trypsin used in step 2 to each well to trypsinize the Matrigel<sup>TM</sup> (125  $\mu$ l trypsin for an 8-well slide, 250  $\mu$ l for a 24-well plate, and 500  $\mu$ l for a 12-well plate).
- 4) Immediately scrape the well to remove the cells and matrigel. For this step, we use a sterile plunger from a 1 ml syringe to scratch the bottom of the well.
- 5) Resuspend the mixture of trypsin, Matrigel<sup>TM</sup>, and cells several times using a P1000 to break up the Matrigel<sup>TM</sup>.
- 6) Collect the trypsin, Matrigel<sup>TM</sup>, and cells in a 15 ml conical tube.
- 7) Repeat steps 3-6 using a fresh volume of trypsin to remove any remaining Matrigel<sup>TM</sup> from the well and collect the mixture in the same 15 ml conical tube.
- 8) Incubate the trypsin, Matrigel<sup>TM</sup>, and cells for 20-30 minutes at 37°C.
- 9) Spin the cells at 150xg in a tissue culture centrifuge for 3 min.
- 10) Resuspend the cell pellet in 5 mls of Resuspension Medium to quench the trypsin (See Media Table for recipe).
- 11) Spin the cells at 150xg in a tissue culture centrifuge for 3 min.
- 12) Resuspend the cell pellet in 1 ml of PBS and transfer to an eppendorf tube.
- 13) Spin the cells for 3 minutes at 2000 rpm in a microcentrifuge. Aspirate the PBS from the cell pellet. (The cell pellet can be flash frozen in liquid nitrogen and stored at –80°C, if necessary. For lysates harvested over several different days, we freeze the pellets on individual days and lyse the cells once all days have been harvested).
- 14) Resuspend pellet in RIPA or NP40 lysis buffer supplemented with protease inhibitors (approximately 100  $\mu$ l/4 wells of a chamber slide).
- 15) Incubate lysate on ice for 15-30 minutes.
- 16) Spin the lysates for 15-20 minutes at 4°C at 14,000xg (maximum speed) in a microcentrifuge to clear the lysates.
- 17) Collect supernatant in new eppendorf tube, and flash freeze the lysate in liquid nitrogen.
- 18) Store at –80°C. Standard BCA or Bradford Assays can be used to quantitate protein levels prior to immunoblotting.