

### ***Harvesting RNA from 3-D Acinar Cultures***

- 1) Aspirate the media from the wells for harvesting.
- 2) Add RNA-STAT-60<sup>TM</sup> (Tel-Test) to each well. For a single well of an 8-well chamber slide, use 250  $\mu$ l of RNA-STAT-60<sup>TM</sup>. For cultures grown in 24-well plates, use 500  $\mu$ l of RNA-STAT-60<sup>TM</sup> per well.
- 3) Immediately scrape the well to remove the cells and matrigel. For this step, we use a sterile plunger from a 1 ml syringe to scratch the bottom of the well.
- 4) Resuspend the mixture of RNA-STAT-60<sup>TM</sup>, Matrigel<sup>TM</sup>, and cells several times using a P1000 to break up the Matrigel<sup>TM</sup>.
- 5) Collect the mixture in an eppendorf tube (approximately 4 wells of a chamber slide per tube).

Following steps are according to manufacturer's instruction:

- 6) Incubate the mixture 5 minutes at room temperature.
- 7) Add 200  $\mu$ l chloroform per 1 ml of RNA-STAT-60<sup>TM</sup> used to the mixture. (Use 200  $\mu$ l for 4 wells of an 8-well chamber slide).
- 8) Shake the tube for 15 seconds in your hand, incubate 2-3 minutes at room temperature.
- 9) Spin the tubes at maximum speed in a microcentrifuge for 15 minutes at 4°C.
- 10) Transfer the upper, aqueous phase to a new eppendorf tube.
- 11) Add 500  $\mu$ l isopropanol per 1 ml RNA-STAT-60<sup>TM</sup> used to the aqueous phase.
- 12) Precipitate the RNA by placing the tubes at 4°C for 30 minutes.
- 13) Spin the tubes at maximum speed in a microcentrifuge for 10 minutes at 4°C to pellet the RNA.
- 14) Wash the pellet with 1 ml 75% ethanol.
- 15) Spin the tubes at 9000xg in a microcentrifuge for 5 minutes.
- 16) Aspirate ethanol and spin once again at 9000xg in a microcentrifuge for 2 minutes to remove excess liquid.
- 17) Dry the pellet at room temperature
- 18) Rehydrate the RNA pellet for 15 minutes at room temperature with approximately 20  $\mu$ l DEPC water (for 4 wells of an 8 well chamber slide).
- 19) Store samples at -80°C.