

***Indirect Immunofluorescent Staining of MCF-10A Acini Cultured In Matrigel
(Detailed)***

- 1) **Fixation:** Aspirate the medium from each well of the chamber slide and immediately fix acini with 2% formalin (Sigma) or 3-4% paraformaldehyde for 20 minutes at room temperature. Acini can also be fixed using a 1:1 mix of methanol:acetone at -20°C for 10-12 minutes. For 8-well chamber slides, 300-400 μl volumes are appropriate for the fixation, permeabilization, and all subsequent washing steps in this protocol. Once fixed, slides can be stored at 4°C for up to two weeks. Rinse wells once in PBS and store in 500 μl /well PBS. Wrap wells in Saran wrap to avoid drying out.
- 2) **Permeabilization:** If formalin or paraformaldehyde is used for fixation, permeabilize with PBS containing 0.5% Triton X-100 for 10 minutes at 4°C . Depending on the antibody utilized for immunostaining, the detergent concentration or duration of permeabilization may require modification.
- 3) **Glycine Rinse:** Rinse 3 times with PBS/Glycine (130 mM NaCl; 7 mM Na_2HPO_4 ; 3.5 mM NaH_2PO_4 ; 100 mM glycine), 10-15 minutes per wash at room temperature.
- 4) **Primary Block:** Incubate with 200 μl /well of IF Buffer (130 mM NaCl; 7 mM Na_2HPO_4 ; 3.5 mM NaH_2PO_4 ; 7.7 mM NaN_3 ; 0.1% BSA; 0.2% Triton X-100; 0.05% Tween-20) +10% goat serum for 45-60 minutes at room temperature.
- 5) **OPTIONAL Secondary Block:** Aspirate the primary block and incubate with 100 μl /well of secondary block (IF Buffer + 10% goat serum + 20 $\mu\text{g}/\text{ml}$ of goat anti-mouse $\text{F}(\text{ab}')_2$ fragment (Jackson ImmunoResearch#115-006-006) for 30-40 minutes.
- 6) **Primary antibody:** Incubate with primary antibody in block solution overnight (15-18 hours) at 4°C . Although optimal antibody concentrations should be determined empirically on a case-by-case basis, a 1:100 to 1:200 dilution of the primary antibody is a good starting point. (See Table 2 for antibodies known to work in this application).
Occasionally, overnight incubation at 4°C elicits liquefaction of the basement membrane and extensive lifting of the acini during subsequent washing steps. This varies with each lot of MatrigelTM used for the morphogenesis assay, and unfortunately, cannot easily be predicted. If this problem arises, it is advisable to perform the primary antibody incubations overnight at room temperature rather than 4°C .
- 7) Rinse 3 times (20 minutes each) with IF Buffer at room temperature with gentle rocking.
- 8) **Secondary Antibody:** Incubate with fluorescent conjugated secondary antibody in IF Buffer + 10% goat serum for 40-50 min at room temperature. We recommend AlexaTM conjugated, highly cross-absorbed secondary antibodies from Molecular Probes used at 1:200 dilution; in our experience, these secondary reagents exhibit low levels of background, and minimal cross-reactivity between species, making them useful for double immunostaining procedures.
- 9) Rinse 1 time (20 minutes) with IF Buffer at room temperature with gentle rocking.
- 10) Rinse 2-3 times with PBS (10 minutes).
- 11) In order to counterstain nuclei, incubate with PBS containing 5 μM TOPRO-3 (Molecular Probes) and/or 0.5 ng /ml 4',6-diamidino-2-phenylindole (DAPI, Sigma) for 15 minutes at room temperature.

- 12) Rinse 1 time with PBS for 5 minutes at room temperature.
- 13) Mount with freshly prepared Prolong Antifade Reagent (Molecular Probes) and allow it dry O/N at RT. Once dry, slides can be stored at 4°C for up to one week or at -20°C for up to two months.